

LACTALIS
INGREDIENTS
PHARMA

Trust our Pharmaceutical lactose



We would be pleased to know
you and your needs.

**Contact us for any specific
need or to receive a sample.**

For more information on our products
www.pharma.lactalisingredients.com
To contact us
lactalisingredientspharma@fr.lactalis.com



Follow us on LinkedIn
Lactalis Ingredients Pharma



LACTALIS
INGREDIENTS
PHARMA

LOW
NITRITE

Lactalpha 200

LACTOSE MONOHYDRATE

200 MESH



Lactalis Ingredients Pharma

Your French trusted pharmaceutical lactose partner

Lactalis Ingredients Pharma contributes to the security of the drug supply chain by delivering the best of lactose to our customers.

Benefiting from the experience and technology of decades of commitment and passion, **we are able today to enter the excipients' market as a trusted pharmaceutical lactose.**

EXPERTISE

QUALITY


RELIABILITY

Lactalpa 200

MILLED LACTOSE MONOHYDRATE



SEM image x 250



LOW
NITRITE
≤ 0.1 ppm

Lactalpa is a range of pharmaceutical grade milled and sieved lactose monohydrate. Lactalpa can easily be used for your OSD applications, in tablets from wet or dry granulation, capsules and sachets. Lactalpa 200 is a 200 mesh pharmaceutical milled lactose conforming with the lactose monohydrate monographs of the current Pharmacopoeias (Ph. Eur., USP-NF and JP).

Production site
Retiers, France

Lactose properties
Diluent / Good mixing
/ Tablet binder

Process
Wet granulation, dry
granulation

Galenic applications
Tablet


Lactalpa 200

Key physical figures

Characteristics	Specification
Bulk density	590 g/L
Tapped density	860 g/L
Hausner ratio	1.46
Carr's index	31%

Particle Size Distribution

PSD by laser diffraction (indicative values)



Laser diffraction

D10: 5 μm
D50: 40 μm
D90: 110 μm


Air jet sieve

<45 μm: 50-65%
<100 μm: ≥90%
<150 μm: ≥96%
<250 μm: ≥98%

— volume
— cumulative volume

Low water absorption for high stability

Sorption isotherm at 25°C



Lactalpa 200 sorption isotherm at 25°C, measured by SPSx-IQ vapor sorption analyzer.

At 25°C, Lactalpa 200 absorbs little water up to 90% relative humidity, and thus shows a great stability.

.....● desorption
.....● sorption